

KPL PhosphaGLO Reserve™ AP Substrate

<u>Catalog No.</u>	<u>Size</u>
5430-0052 (55-60-01)	30 mL
5430-0053 (55-60-02)	100 mL

DESCRIPTION

KPL PhosphaGLO Reserve™ AP Substrate contains a dioxetane-based chemiluminescent substrate designed for use with phosphatase-labeled (AP) reporter molecules. KPL PhosphaGLO Reserve offers great improvements in the way of signal intensity. These products are specifically designed for the detection of proteins that are (1) difficult to detect because they are in such low quantities or (2) are from samples that are precious.

KPL PhosphaGLO Reserve AP Substrate is provided as a stable one-component solution. The product provides rapid and accurate identification of proteins that are of low abundance and potentially limited availability. Given the increased sensitivity, less target may be required on Western blots and ELISAs.

Two sizes are available. Results can be obtained on X-ray film or a chemiluminescent imager to provide a permanent record. A luminometer should be used for ELISA.

NOTE: This product may arrive frozen.

CONTENTS

5430-0052 (55-60-01) contains:
1 x 30 mL KPL PhosphaGLO Reserve AP Substrate

5430-0053 (55-60-02) contains:
1 x 100 mL KPL PhosphaGLO Reserve AP Substrate

STORAGE/STABILITY

Store this product at 2–8°C. If it arrives frozen, allow it to thaw prior to use. KPL PhosphaGLO Reserve Solution should remain stored in its original container and protected from light.

The product is stable for a minimum of two years from date of receipt when stored under proper conditions.

KPL PHOSPHAGLO RESERVE CHEMILUMINESCENT SUBSTRATE USER'S GUIDE

- ⇒ KPL PhosphaGLO Reserve can be used with nitrocellulose and PVDF membranes. For best results, nitrocellulose is recommended.
- ⇒ The KPL PhosphaGLO Reserve working solution should be protected from light and warmed to room temperature prior to use.
- ⇒ KPL PhosphaGLO Reserve is an extremely sensitive substrate. Insufficient washing of membranes or contamination of substrate with AP will result in non-specific background.
- ⇒ Because of KPL PhosphaGLO Reserve's super sensitivity, the AP conjugate should be titrated to give the best results.
- ⇒ Do not allow KPL PhosphaGLO Reserve to contact the film. If this occurs, KPL PhosphaGLO Reserve solution will cause dark spots to appear on the film.
- ⇒ KPL PhosphaGLO Reserve emits light over the course of 5-7 days. Because of its high light intensity, images may be captured over the course of many different intervals. Initially, perform a 1 minute and a 10 minute exposure.

APPLICATIONS

KPL PhosphaGLO Reserve has been optimized for Western blotting and dot blotting applications. It is also suitable for use in microwell applications such as ELISA. The following is a recommended procedure for Western blot detection.

WESTERN BLOT DETECTION

Suggested Reagents/Equipment Not Included

1. Primary antibody
2. AP-labeled secondary
3. Nitrocellulose or PVDF membrane
4. Blocking Solution (See RELATED PRODUCTS)
5. X-ray film (double emulsion) or CCD Imager
6. Platform shaker or rocker
7. Developing chemicals/equipment
8. Incubation trays or tubes
9. Wash solution (see RELATED PRODUCTS)
10. Assay buffer (see RELATED PRODUCTS)

KPL PhosphaGLO Reserve™ AP Substrate


<u>Catalog No.</u>	<u>Size</u>
5430-0052 (55-60-01)	30 mL
5430-0053 (55-60-02)	100 mL

CONJUGATE OPTIMIZATION PRIOR TO DETECTION

- ⇒ Before beginning the assay, the optimal conjugate dilution should be determined. The use of highly sensitive chemiluminescent substrates on Western blots can cause high background if the conjugate concentration is not optimized. Each lot of conjugate should be optimized, as slight differences in activity can produce major differences in background.
- ⇒ Recommended conjugate dilutions should be tested at a range from 1/10,000 to 1/100,000 of a 0.1 mg/mL stock.

WESTERN BLOT DETECTION AT A GLANCE

Total time: 4 hours



KPL PhosphaGLO Reserve™ AP Substrate

Catalog No.
 5430-0052 (55-60-01)
 5430-0053 (55-60-02)

Size
 30 mL
 100 mL

STEPS	CRITICAL POINTS
1. Block the membrane by immersing in block solution (1X Detector Block is recommended) using a minimum of 0.2 mL/cm ² of membrane. Block at room temperature for 1 hour with gentle rocking or shaking, or stationary at 2–8°C overnight.	<i>Example: for a 10 x 10 cm blot, use 20 mL of block solution. Make sure to use a container of proper size that allows the block solution to float freely over the membrane.</i>
2. Incubate membrane with primary antibody or serum sample for at least 1 hour. This antibody should be added directly to the Block Solution that was used for blocking (Step 1).	<i>Test serial dilutions through a dot blot to determine the optimal working dilution. Alternatively, use the recommended concentration determined by the primary antibody supplier.</i>
3. Wash the membrane in a generous amount of 1X Wash Solution (at least 25 mL for a 100 cm ² membrane). Wash membrane 3 times for 5 minutes each, followed by one 10-minute wash.	<i>1X TBS-TWEEN may also be used.</i>
4. Dilute appropriate conjugate 1/10,000 – 1/100,000 (of a 0.1 mg/mL stock) in freshly prepared conjugate diluent using a minimum of 0.2 mL/cm ² of membrane.	<i>Example: 2 µL conjugate + 20 mL diluent. Suggested diluents include KPL Detector Block and TBS-TWEEN. The optimal dilution may vary for different lots of conjugate.</i>
5. Incubate blot with diluted conjugate for one hour at room temperature.	<i>Titrate the conjugate to determine the optimal working dilution.</i>
6. During the conjugate incubation step, prepare 0.05 mL/cm ² of membrane to be detected.	<i>Prepare KPL PhosphaGLO Reserve in advance to allow it to come to room temperature prior to its use. Cover with foil to minimize light exposure.</i>

STEPS	CRITICAL POINTS
7. After the conjugate incubation, wash as described in step 3.	
8. Pour off the remaining wash buffer from the blot and place the membrane on a sheet protector or a dry tray.	
9. Rinse membrane 2 x 2 minutes with 1X Assay buffer or 0.2 M Tris, pH 9.5-9.7.	
10. Gently pipette 0.05 mL/cm ² of previously prepared KPL PhosphaGLO Reserve over the entire membrane. Incubate without rocking for 1 minute.	<i>Example: for a 10 x 10 cm blot, use 5 mL of KPL PhosphaGLO Reserve. The surface tension of the substrate will keep it on the surface of the membrane.</i>
11. Lift the membrane with forceps and blot the excess substrate onto a piece of filter paper. Seal the membrane in clear plastic and expose to X-ray film for 10 seconds to 1 minute. Adjust exposure time as needed.	<i>Excessive substrate on the blot will contribute to background. Take caution to ensure the surface of the membrane to which the assay reagents were applied is facing the film. Do not allow the film to get wet, nor move during exposure.</i>
12. <i>Optional:</i> Chemiluminescent Imager Detection. Incubate the blot for twice the time typically used for film. If the imager provides stacking capabilities, capture exposures at 5 minute intervals for 1 hour to maximize signal. The optimal exposure can be chosen.	<i>Optimal exposure time should be determined by the signal to noise ratio and the amount of conjugate used. When using greater amounts of conjugate, 10 seconds may provide acceptable results.</i>
	<i>Follow the manufacturer's recommendations regarding the set up and operation of the imager.</i>

KPL PhosphaGLO Reserve™ AP Substrate

Catalog No.
 5430-0052 (55-60-01)
 5430-0053 (55-60-02)

Size
 30 mL
 100 mL

TROUBLESHOOTING

Problem 1: No Signal

Possible Cause	Corrective Measure
<ul style="list-style-type: none"> Inactive alkaline phosphatase 	Verify enzyme activity by mixing 10 µL of diluted conjugate with 1 mL of substrate (the substrate should glow in the dark).
<ul style="list-style-type: none"> No binding of conjugate to the primary antibody 	Confirm correct specificity of the conjugate for the primary antibody; i.e. no AP-anti-rabbit with a mouse primary antibody.
<ul style="list-style-type: none"> No transfer of target to membrane 	Use a protein stain on unblocked membrane to verify attachment of target protein or use a pre-stained protein marker to monitor transfer.
<ul style="list-style-type: none"> Detection of non-blotted side of membrane 	Ensure correct orientation of the membrane during the assay and film exposure.
<ul style="list-style-type: none"> Missed step in procedure 	Review procedure to ensure all steps were followed.

Problem 2: Weak Signal

Possible Cause	Corrective Measure
<ul style="list-style-type: none"> Insufficient amount of antibody 	Optimize antibody concentrations. Affinity of the primary antibody may change after proteins are denatured through SDS-PAGE.
<ul style="list-style-type: none"> Insufficient protein loaded or transferred 	Increase the amount of protein loaded onto the gel.
<ul style="list-style-type: none"> Insufficient incubation of primary antibody 	Increase the incubation times for weak primary antibodies.
<ul style="list-style-type: none"> Insufficient exposure time 	Increase the time of exposure to film.
<ul style="list-style-type: none"> Excessive washing beyond recommended procedure 	Follow the procedure as written.

Problem 3: Excessive signal, nonspecific bands or general background

Possible Cause	Corrective Measure
<ul style="list-style-type: none"> Overexposure of film to signal 	Expose the membrane to film for a shorter period of time.
<ul style="list-style-type: none"> Insufficient blocking or washing 	Increase blocking and washing time or increase number of washes. Vary type of block used.
<ul style="list-style-type: none"> Excessive antibody used for detection 	Optimize conjugate concentration. Reduce antibody concentrations; optimal conjugate dilution should be 1/10,000 – 1/100,000 of a 0.1 mg/mL stock. Or, decrease the amount of primary antibody.
<ul style="list-style-type: none"> Excessive protein loaded on the gel 	Decrease the amount of protein loaded onto the gel.

Problem 4: Poorly Defined or “Fuzzy” Bands or Dots

Possible Cause	Corrective Measure
<ul style="list-style-type: none"> Poor transfer of protein to membrane 	Follow manufacturer’s recommended procedure or contact the manufacturer for additional support regarding the blotting apparatus.
<ul style="list-style-type: none"> Excessive substrate 	Remove excess substrate before exposure of the membrane to film.
<ul style="list-style-type: none"> Ghost images from shifted position of film during development 	Avoid movement of film over membrane during exposure period.
<ul style="list-style-type: none"> Inadequate handling of membranes 	Certain membranes require special handling. Check with the membrane vendor for correct procedures.

KPL PhosphaGLO Reserve™ AP Substrate

<u>Catalog No.</u>	<u>Size</u>
5430-0052 (55-60-01)	30 mL
5430-0053 (55-60-02)	100 mL

Stripping and Re-probing a Western Blot

This protocol is adapted from Kaufmann, *et. al.*¹¹. After performing protein transfer, detection with KPL PhosphaGLO Reserve and film exposure, membranes may be stripped and re-probed with new primary and secondary antibodies.

1. Strip antibodies by incubating blot for 30 - 90 minutes at 70°C in erasure buffer: 2% SDS (w/v), 62.5 mM Tris-HCl (pH 6.8 at 20°C), 100 mM β-mercaptoethanol.
2. Wash 2 times, for 10 minutes each, in TBS: 10 mM Tris-HCl (pH 7.4 at 20°C), 150 mM NaCl.
3. Block for 2.5 hours in Block Solution.
4. Repeat detection procedure.

PRODUCT SAFETY AND HANDLING

See SDS (Safety Data Sheet) for this product.

RELATED PRODUCTS	CAT. NO.
KPL PhosphaGLO™ AP Substrate	5430-0054 (55-60-03)
KPL Detector™ Block	5920-0004 (71-83-00)
KPL Milk Diluent/Blocking Solution	5140-0011 (50-82-01)
KPL 10% BSA Diluent/Blocking Solution	5140-0006 (50-61-00)
KPL Wash Solution	5150-0008 (50-63-00)
KPL Coating Solution	5150-0014 (50-84-00)
KPL Phosphatase Assay Buffer	5960-0017 (50-63-14)

PhosphaGLO and PhosphaGLO Reserve are trademarks of SeraCare Life Sciences, Inc.
 Detector is a trademark of SeraCare Life Sciences, Inc.
 TWEEN is a trademark of ICI Americas, Inc.

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.

REFERENCES

1. Kricka, L. J. (1991). Chemiluminescent and Bioluminescent Techniques. *Clin. Chem.* 37(9): 1472 - 1481.
2. Knecht, D.A. and R.L. Dimond (1984). Visualization of Antigenic Proteins on Western Blots. *Anal. Biochem.* 136: 180 - 184.
3. Blake, M.S., et al (1984). A Rapid, Sensitive Method For Detection of Alkaline Phosphatase Conjugated Antibody on Western Blots. *Anal. Biochem.* 136: 175 - 178.
4. Isacsson, V. and G. Wettermark (1974). Chemiluminescence in Analytical Chemistry. *Anal. Chim. Acta.* 68: 339 - 362.
5. Towbin, H., T. Staehelin and J. Gordon (1979). Electrophoretic Transfer of Proteins From Polyacrylamide Gels to Nitrocellulose Sheets: Procedure and Some Applications. *Proc. Natl. Acad. Sci. USA.* 76: 4350 - 4354.
6. Bittner, M., P. Kupferer and C. F. Morris (1980). Electrophoretic Transfer of Proteins and Nucleic Acids From Slab Gels to Diazobenzoyloxymethyl Cellulose or Nitrocellulose Sheets. *Anal. Biochem.* 102: 459 - 471.
7. Burnette, W.N. (1980). "Western Blotting": Electrophoretic Transfer of Proteins From Sodium Dodecyl Sulfate-Polyacrylamide Gels to Unmodified Nitrocellulose or Nitrocellulose Sheets. *Anal. Biochem.* 112: 195 - 203.
8. Reinhart, M.P. and D. Malamud (1982). Protein Transfer From Isoelectric Focusing Gels: The Native Blot. *Anal. Biochem.* 123: 229 - 235.
9. Gooderham, K (1983). Protein Blotting. In J. Walker and W. Gaastra (eds.), *Techniques in Molecular Biology*. Croom Helm Ltd. Publishers, London.
10. Southern, E.M. (1975). Detection of Specific Sequences Among DNA Fragments Separated by Gel Electrophoresis. *J. Mol. Biol.* 98: 503 - 517.
11. Kaufmann, Ewing and Shaper (1987). The Erasable Western Blot. *Anal. Biochem.* 161, 89 - 95.